

Pipettor Compatibility Chart

FILTERED, STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	20 µL	100 µL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL
					4117NSFS	4137NSFS	4237NAFS	4237NBFS	4237NSFS	4287NSFS	4927NSFS	4337NSFS	4347NSFS
Brand	Transferpette S	1 µL	0.1 - 1 µL	SINGLE	4117NSFS	4137NSFS							
		2.5 µL	0.1 - 2.5 µL	SINGLE	4117NSFS	4137NSFS							
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS							
		20 µL	2 - 20 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS				
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS				
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
		1000 µL	100 - 1000 µL	SINGLE								4337NSFS	4347NSFS
		200 µL	20 - 200 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
		300 µL	15 - 300 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
Capp	Aero	2 µL	0.2 - 2 µL	SINGLE	4117NSFS	4137NSFS	4237NAFS	4237NBFS					
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS							
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS		4927NSFS		
		50 µL	5 - 50 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
		200 µL	25 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
		300 µL	50 - 300 µL	SINGLE							4927NSFS		
		1000 µL	100 - 1000 µL	SINGLE								4337NSFS	4347NSFS
		10 µL	0.5 - 10 µL	MULTI	4117NSFS	4137NSFS							
	50 µL	5 - 50 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
	Bravo	2 µL	0.1-2 µL	SINGLE	4117NSFS	4137NSFS	4237NAFS	4237NBFS					
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS							
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS		4927NSFS		
		50 µL	5 - 50 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
		1000 µL	100 - 1000 µL	SINGLE								4337NSFS	4347NSFS
		10 µL	0.5 - 10 µL	MULTI	4117NSFS	4137NSFS							
50 µL		5 - 50 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
Eppendorf	Reference	2.5 µL	0.1 - 2.5 µL	SINGLE	4117NSFS	4137NSFS							
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS							
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS				
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS				
		1000 µL	100 - 1000 µL	SINGLE								4337NSFS	4347NSFS
	Reference 2	2.5 µL	0.1 - 2.5 µL	SINGLE	4117NSFS	4137NSFS							
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS							
		20 µL	2 - 20 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS		4927NSFS		

Pipettor Compatibility Chart

FILTERED, STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	20 µL	100 µL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117NSFS	4137NSFS	4237NAFS	4237NBFS	4237NSFS	4287NSFS	4927NSFS	4337NSFS	4347NSFS	
Eppendorf	Reference 2	200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		300 µL	30 - 300 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE								4337NSFS	4347NSFS	
	Research	2.5 µL	0.1 - 2.5 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS							
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS					
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		300 µL	30 - 300 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
		10 µL	0.5 - 10 µL	MULTI	4117NSFS	4137NSFS								
		100 µL	10 - 100 µL	MULTI			4237NAFS	4237NBFS	4237NAFS					
	300 µL	30 - 300 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS				
	Research plus	2.5 µL	0.1 - 2.5 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL (yellow)	SINGLE			4237NAFS		4237NAFS					
		20 µL	2 - 20 µL (light gray)	SINGLE	4117NSFS	4137NSFS								
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS					
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		300 µL	30 - 300 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
		10 µL	0.5 - 10 µL	MULTI	4117NSFS	4137NSFS								
		300 µL	20 - 300 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
	Xplorer / Xplorer Plus	10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	1 - 20 µL	SINGLE			4237NAFS		4237NAFS	4287NSFS				
		100 µL	5 - 100 µL	SINGLE				4237NBFS	4237NAFS	4287NSFS				
		200 µL	10 - 200 µL	SINGLE					4237NAFS	4287NSFS				
		300 µL	15 - 300 µL	SINGLE							4927NSFS			
1000 µL		50 - 1000 µL	SINGLE									4337NSFS	4347NSFS	
10 µL		0.5 - 10 µL	MULTI	4117NSFS	4137NSFS									
100 µL		5 - 100 µL	MULTI				4237NBFS	4237NAFS	4287NSFS					
300 µL		15 - 300 µL	MULTI							4927NSFS				
1200 µL		50 - 1200 µL	MULTI										4347NSFS	
Finnpipette		4500	10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS							
	20 µL		2 - 20 µL	SINGLE			4237NAFS		4237NAFS		4927NSFS			
	50 µL		5 - 50 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
	200 µL		20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			

Pipettor Compatibility Chart

FILTERED, STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	20 µL	100 µL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117NSFS	4137NSFS	4237NAFS	4237NBFS	4237NSFS	4287NSFS	4927NSFS	4337NSFS	4347NSFS	
Finnpipette	4500	1000 µL	200 - 1000 µL	SINGLE								4337NSFS	4347NSFS	
		50 µL	5 - 50 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
	F1	2 µL	0.2 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		5 µL	0.5 - 5 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL (pink)	SINGLE	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL (yellow)	SINGLE			4237NAFS		4237NAFS		4927NSFS			
		20 µL	2 - 20 µL (micro)	SINGLE	4117NSFS	4137NSFS								
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
	300 µL	30 - 300 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS				
	F2	2 µL	0.2 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		5 µL	0.5 - 5 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS		4927NSFS			
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
		100 µL	10 - 100 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		300 µL	30 - 300 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
	Focus	50 µL	5 - 50 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS		4927NSFS			
		300 µL	30 - 300 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
	Novus	10 µL	1 - 10 µL (pink)	SINGLE	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL (yellow)	SINGLE			4237NAFS		4237NAFS		4927NSFS			
		50 µL	5 - 50 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS		4927NSFS			
		300 µL	30 - 300 µL	SINGLE					4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
10 µL		1 - 10 µL	MULTI	4117NSFS	4137NSFS									
50 µL		5 - 50 µL (yellow)	MULTI			4237NAFS	4237NBFS	4237NAFS		4927NSFS				
300 µL		30 - 300 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS				
1200 µL	100 - 1200 µL	MULTI									4337NSFS	4347NSFS		
Gilson	Pipetman Classic	2 µL	0.2 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS					
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS

Pipettor Compatibility Chart

FILTERED, STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	20 µL	100 µL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117NSFS	4137NSFS	4237NAFS	4237NBFS	4237NSFS	4287NSFS	4927NSFS	4337NSFS	4347NSFS	
Gilson	Pipetman Classic	20 µL	2 - 20 µL	MULTI			4237NAFS		4237NAFS	4287NSFS	4927NSFS			
		200 µL	20 - 200 µL	MULTI				4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		300 µL	30 - 300 µL	MULTI					4237NAFS	4287NSFS	4927NSFS			
	Pipetman L/M/G	2 µL	0.2 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS			4927NSFS		
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS			4927NSFS		
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
		10 µL	0.5 - 10 µL	MULTI	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	MULTI			4237NAFS		4237NAFS			4927NSFS		
		200 µL	20 - 200 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
	300 µL	20 - 300 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS				
	Pipetman Neo	2 µL	0.2 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS					
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE					4237NAFS	4287NSFS	4927NSFS			
	Pipetman Ultra	10 µL	1 - 10 µL	SINGLE		4137NSFS								
20 µL		2 - 20 µL	SINGLE			4237NAFS		4237NAFS			4927NSFS			
100 µL		20 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS			4927NSFS			
300 µL		20 - 300 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS				
1000 µL		100 - 1000 µL	SINGLE									4337NSFS	4347NSFS	
Hamilton	Softgrip	2 µL	0.2 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL	SINGLE		4137NSFS								
		25 µL	2.5 - 25 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		300 µL	30 - 300 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
		300 µL	30 - 300 µL	MULTI				4237NBFS	4237NAFS	4287NSFS	4927NSFS			
Hirschmann	labopette	3 µL	0.1 - 3 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS			4927NSFS		
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS			4927NSFS		
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS

Pipettor Compatibility Chart

FILTERED, STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	20 µL	100 µL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117NSFS	4137NSFS	4237NAFS	4237NBFS	4237NSFS	4287NSFS	4927NSFS	4337NSFS	4347NSFS	
HTL Pipette Discovery	Pro / Comfort / LABMATE pro / Optipette	2 µL	0.2 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS			4927NSFS		
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
IKA PETTE	vario	2 µL	0.1 - 2.0 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS					
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS					
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE								4337NSFS	4347NSFS	
		10 µL	0.5 - 10 µL	MULTI	4117NSFS	4137NSFS								
		100 µL	10 - 100 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
Nichiryo	Oxford Benchmate / Benchmate II	2 µL	0.1 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS		4927NSFS			
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE								4337NSFS	4347NSFS	
	EX II MULTI	10 µL	0.5 - 10 µL	MULTI	4117NSFS	4137NSFS								
		100 µL	5 - 100 µL	MULTI			4237NAFS	4237NBFS	4237NAFS		4927NSFS			
		200 µL	20 - 200 µL	MULTI				4237NBFS	4237NAFS		4927NSFS			
		300 µL	30 - 300 µL	MULTI						4287NSFS	4927NSFS			
	Nichipet EX / EX II / EX PLUS II / LE / Premium / Premium LT	10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS		4927NSFS			
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
1000 µL		100 - 1000 µL	SINGLE								4337NSFS	4347NSFS		
Rainin	Pipet-Lite XLS (Universal)	2 µL	0.1 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20µL	SINGLE			4237NAFS		4237NAFS		4237NSFS			
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS		4237NSFS			
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4237NBFS	4237NSFS			
		300 µL	20 - 300 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4237NBFS	4237NSFS			
		1000 µL	100 - 1000 µL	SINGLE								4927NSFS	4287NSFS	

Pipettor Compatibility Chart

FILTERED, STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	20 µL	100 µL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117NSFS	4137NSFS	4237NAFS	4237NBFS	4237NSFS	4287NSFS	4927NSFS	4337NSFS	4347NSFS	
Sartorius	(Biohit) mLINE	3 µL	0.1 - 3 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.5-10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS		4927NSFS			
		100 µL	5 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 -1000 µL	SINGLE									4337NSFS	4347NSFS
		10 µL	0.5 -10 µL	MULTI	4117NSFS	4137NSFS								
		100 µL	5 - 100 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		300 µL	30 - 300 µL	MULTI				4237NBFS	4237NAFS	4287NSFS	4927NSFS			
	(Biohit) Proline Plus / Proline	3 µL	0.1 - 3 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS		4927NSFS			
		50 µL	5 - 50 µL	SINGLE			4237NAFS		4237NAFS					
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS				
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
		10 µL	0.5 -10 µL	MULTI	4117NSFS	4137NSFS								
		100 µL	10 - 100 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS				
	300 µL	30 - 300 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS				
	Picus	2 µL	0.2 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.2 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		120 µL	5 - 120 µL	SINGLE						4287NSFS	4927NSFS			
		300 µL	10 - 300 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	50 - 1000 µL	SINGLE									4337NSFS	4347NSFS
		10 µL	0.2 - 10 µL	MULTI	4117NSFS	4137NSFS								
		120 µL	5 - 120 µL	MULTI				4237NBFS		4287NSFS	4927NSFS			
		300 µL	10 -300 µL	MULTI	4117NSFS	4137NSFS	4237NAFS	4237NBFS						
300 µL		10 - 300 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS				
1200 µL		50 - 1200 µL	MULTI									4337NSFS	4347NSFS	
Tacta	3 µL	0.1 - 3 µL	SINGLE	4117NSFS	4137NSFS									
	10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS									
	20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS	4287NSFS					
	100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS					
	200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS					
	1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS	
	3 µL	0.1 - 3 µL	MULTI	4117NSFS	4137NSFS									
10 µL	0.5 - 10 µL	MULTI	4117NSFS	4137NSFS										

Pipettor Compatibility Chart

FILTERED, STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	20 µL	100 µL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117NSFS	4137NSFS	4237NAFS	4237NBFS	4237NSFS	4287NSFS	4927NSFS	4337NSFS	4347NSFS	
Sartorius	Tacta	100 µL	5 - 100 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS				
		300 µL	30 - 300 µL	MULTI				4237NBFS	4237NAFS	4287NSFS				
Socorex	Acura	2 µL	0.1 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS					
		50 µL	5 - 50 µL	SINGLE			4237NAFS		4237NAFS					
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS		4927NSFS			
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
		10 µL	0.5 - 10 µL	MULTI	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL	MULTI	4117NSFS	4137NSFS								
		50 µL	5 - 50 µL	MULTI			4237NAFS		4237NAFS					
		100 µL	10 - 100 µL	MULTI			4237NAFS	4237NBFS	4237NAFS		4927NSFS			
		200 µL	20 - 200 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS				
	350 µL	40 - 350 µL	MULTI				4237NBFS	4237NAFS	4287NSFS	4927NSFS				
	Calibra 822	2 µL	0.2 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	1 - 10 µL (200µL tips)	SINGLE			4237NAFS		4237NAFS					
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS				
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS				
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
10 µL		1 - 10 µL	MULTI	4117NSFS	4137NSFS									
100 µL		10 - 100 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS					
200 µL	20 - 200 µL	MULTI				4237NBFS	4237NAFS	4287NSFS						
VWR	Ergonomic	2 µL	0.2 - 2 µL	SINGLE	4117NSFS	4137NSFS								
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS								
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS					
		50 µL	5 - 50 µL	SINGLE			4237NAFS		4237NAFS	4287NSFS				
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS				
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS				
		250 µL	50 - 250 µL	SINGLE				4237NBFS	4237NAFS	4287NSFS	4927NSFS			
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS	4347NSFS
		10 µL	1 - 10 µL	MULTI	4117NSFS	4137NSFS								
		50 µL	5 - 50 µL	MULTI			4237NAFS		4237NAFS	4287NSFS				
		200 µL	20 - 200 µL	MULTI			4237NAFS	4237NBFS	4237NAFS	4287NSFS				
		300 µL	50 - 300 µL	MULTI				4237NBFS	4237NAFS	4287NSFS	4927NSFS			

Pipettor Compatibility Chart

FILTERED, STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	20 µL	100 µL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL
					4117NSFS	4137NSFS	4237NAFS	4237NBFS	4237NSFS	4287NSFS	4927NSFS	4337NSFS	4347NSFS
VWR	Ultra	2 µL	0.1 - 2 µL	SINGLE	4117NSFS	4137NSFS							
		10 µL	0.5 - 10 µL	SINGLE	4117NSFS	4137NSFS							
		20 µL	2 - 20 µL	SINGLE			4237NAFS		4237NAFS		4927NSFS		
		100 µL	10 - 100 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
		200 µL	20 - 200 µL	SINGLE			4237NAFS	4237NBFS	4237NAFS	4287NSFS	4927NSFS		
		1000 µL	100 - 1000 µL	SINGLE									4337NSFS